

GELSCANによるSYPRO Ruby 検量線  
 300ppi、ガンマ1.0、16bitTIFF、Redチャンネル  
 file:20071204\_154226.tif

バンド	50kDaバンド タンパク質 量[ng/band]	ImageJで 計測した面 積
1	3000	6984
2	1500	5102
3	750	3536
4	375	2481
5	189	1667
6	93	988
7	48	622
8	23.4	318
9	11.7	232
10	6	157
11	3	86
12	1.5	28
14	1500	4949



